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Vente exclusive en pharmacie

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BP 70103 - 54183 HEILLECOURT - FRANCE

Fabriqué par:  Microlife Corporation

9F, 431, RuiGuang Road, NeiHu

Taipei 11492, Taiwan, R.O.C.

Mandataire: Microlife AG

Espenstrasse 139

9443 Widnau - Suisse

Made in China

CÉ0044
Microlife Corporation Taipei - TW

D I G I T

Digitensio PRO

BP 3BM1-4D

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Automatic Blood Pressure Monitor

Instruction Manual

This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

Read the instructions carefully before using this device.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

Keep dry.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Type BF applied part.

Table of Contents

1. Introduction

- 1.1. Features of BP 3BM1-4D
- 1.2. Important information about self-measurement

2. Important Information on the Subject of Blood Pressure and its Measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4. MAM technology (Average mode)
- 2.5. IHD - Pulse Arrhythmia Detection

3. Components of your Blood Pressure Monitor

4. Using your Blood Pressure Monitor

- 4.1. Inserting the batteries
- 4.2. Using an mains adapter (special accessory)
- 4.3. Cuff connection
- 4.4. Setting the time and date
- 4.5. Select the user
- 4.6. Select Measuring Mode: Standard or Averaging Mode

5. Carrying out a Measurement

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Memory storage and recall of the measurements
- 5.6. Discontinuing a measurement
- 5.7. Delete memories

6. Error Messages/Malfunions

7. Care and Maintenance, Recalibration

8. Guarantee

9. Technical Specifications and reference to standards

Intended use:

This oscillometric blood pressure monitor is intended for measuring non-invasive blood pressure in people aged 12 years or older.

1. Introduction

1.1. Features of your BP 3BM1-4D

Your blood pressure monitor BP 3BM1-4D is a fully automatic, digital blood pressure measuring device for use on the upper arm with integrated MAM and IHD technology. This monitor enables very fast and reliable measurement of the systolic and diastolic blood pressure as well your pulse by use of the oscillometric method.

This device is easy to use and is proven in clinical studies to provide excellent accuracy. The large display makes it easy to visualize operational status and cuff pressure during a measurement.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- The measurement results given by this device is not a diagnosis. It is not replacing the need for the consultation of a physician, especially if not matching the patient's symptoms. Do not rely on the measurement result only, always consider other potentially occurring symptoms and the patient's feedback. Calling a doctor or an ambulance is advised if needed.
- Self-measurement means **Control**, not diagnosis or treatment. Your values must always be discussed with your doctor. **Never alter the dosages of any medication without direction from your doctor.**
- The pulse reading is **not** suitable for checking heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference

Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.

2. Important Information on the Subject of Blood Pressure and its Measurement

2.1. How does high/low blood pressure arise?

Your blood pressure level is determined in the circulatory center of your brain. Your nervous system allows your body to adapt or alter blood pressure in response to different situations. Your body alters your pulse and the width of blood vessels through changes in muscles in the walls of blood vessels.

Your blood pressure reading is highest when the heart pumps or ejects blood – this is called your Systolic Blood Pressure. Your blood pressure reading is lowest when the heart rests (in-between beats) – this is called your Diastolic Blood Pressure.

Note: Maintain blood pressure values within a «normal» range(s) in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood pressure is over 140 mmHg. If you obtain readings in this range, consult your doctor immediately. High blood pressure values over time damage blood vessels, vital organs such as the kidney and even your heart.

When blood pressure values are too low, i.e. systolic values under 100 mmHg and/or diastolic values under 60 mmHg, consult your doctor.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to independently alter drug doses prescribed by your doctor.**

Table for classifying blood pressure values (units mmHg): according to World Health Organization:

Colour	Range	Systolic	Diastolic	Recommendation
Light green	blood pressure too low	↓ 100	↓ 60	Consult your doctor
Light green	blood pressure optimum	100 - 120	60 - 80	self-check
Green	blood pressure normal	120 - 130	80 - 85	self-check
Yellow	blood pressure slightly high	130 - 140	85 - 90	Consult your doctor
Orange	blood pressure too high	140 - 160	90 - 100	Seek medical advice
Red	blood pressure far too high	160 - 180	100 - 110	Seek medical advice
Dark red	blood pressure dangerously high	180 ↑	110 ↑	Urgently seek medical advice!

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Further information

- If your values are mostly «normal» under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». In any case, please discuss the values with your doctor.
- Correctly measured diastolic blood pressure values above 120 mmHg require **immediate medical treatment**.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) over time are associated with considerable risks to health. Blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). This can result in a deficient supply of blood to important organs (heart, brain, muscles). Additionally, the heart will become structurally damaged.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can cause organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are lifestyle changes you can make to prevent and reduce high blood pressure. These measures must be part of a healthy lifestyle and include:

A) Eating habits

- Strive for a normal weight as prescribed by your doctor. Reduce overweight!
- Avoid excessive consumption of common salt. (Many «packaged foods» contain high levels of salt).
- Avoid fatty foods. (Packaged foods are frequently high in fats).

B) Previous illnesses

Consistently follow any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

C) Habits

- Eliminate smoking
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee, tea, chocolate, etc.)

D) Physical conditioning

- After a preliminary medical examination, exercise regularly.
- Choose sports which require endurance and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise activities. Your physician will help you develop an exercise routine that is appropriate for you.

2.4. MAM technology (Average Mode)

- «MAM technology» is a new type of concept for optimum reliability in self-measurement of blood pressure.
- An advanced measurement accuracy is achieved by the automatic analysis of three successive measurements.
- The new system provides reliable values for the doctor and can be used as the basis for reliable diagnostics and medication therapy for high blood pressure.

A) Why MAM?

- Human blood pressure is not stable
- Scattering of devices

B) Key advantages

Reduction of:

- Device scattering
- Insufficient rest prior to measurement
- Movement artefacts
- Cuff positioning influences

C) Medical benefits

- Improved accuracy
- Reliable patient self-measurement data for the doctor
- Safe hypertension diagnostic
- Reliable therapy control

D) Measurement sequence

- With the full measurement cycles, the total measurement time remains less than 3 minutes, compared to a single measurement time of 1.5 minute.
- Single results are not displayed.
- Due to the «Data Analysis» result, a 4th measurement may be applied.

1. measurement

15 seconds resting time

2. measurement

15 seconds resting time

3. measurement

Analysis

Result

2.5. IHD - Pulse Arrhythmia Detection

Appearance of the Arrhythmia indicator

The appearance of this symbol  signifies that a certain pulse irregularity was detected during the measurement. The result can vary from your normal blood pressure. As a rule this is not a cause for concern; however, if the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend you inform your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Arrhythmia indicator

This BP 3BM1-4D device is an oscillometric blood pressure measuring device that also analyses the pulse frequency during measurement. The accuracy of this device has been clinically validated.

If pulse irregularities occur during measurement, the arrhythmia symbol is displayed after the measurement. If the symbol appears **more frequently** (e.g. several times per week on measurements performed daily) **or if it suddenly appears more often than usual**, we recommend the patient to seek medical advice.



The device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

3. Components of the Blood Pressure Monitor

The illustration shows the blood pressure monitor, consisting of:


a) Measuring unit

Single Mode /
MAM – Average Mode

b) Upper arm cuff



M-L cuff for arm circumference 22-42 cm


4. Using your Blood Pressure Monitor

4.1. Inserting the batteries

Insert batteries immediately after unpacking the device. The battery compartment is located on the back of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (4 x size AA 1,5V), observing the indicated polarity.
- If a battery warning  appears in the display, the batteries are almost flat and must be replaced; After battery warning  appears, the device will not work until batteries have been replaced.


Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5 V Batteries. The use of 1.2 V Accumulators is not recommended.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.

4.2. Using an mains adaptor (special accessory)

This blood pressure monitor can be operated with the mains adaptor (output 6 V DC / 600 mA, DIN plug).

- Plug the plug into the socket at the back of the instrument.
- Plug the mains adaptor into a 230 V or 110 V power socket. Test that power is available by pressing the  button.


Note:

- No power is taken from the batteries while the mains adaptor is connected to the instrument.
- If the mains voltage is interrupted during the measurement (e.g. by accidental removal of the mains adaptor from the mains socket), the instrument must be reset by removing the plug from its socket and re-inserting the connections.
- Please consult a specialist dealer if you have questions relating to the mains adaptor.

4.3. Cuff connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.


4.4. Setting the time and date

This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement.


After new batteries have been inserted, the device will request a user to set up time/date.

(Example: Entering 2008-08-15 Time 09:30)

1. After new batteries have been inserted, the display comes into direct year setting, during which the 4 characters (2005) blink.

2. Enter the correct year by pressing «**MEMORY**» button. Each click will change the clock one year. Example: 3 x press.

3. Press «**TIME**» button. The display now switches to the current date, during which the first character (month) blinks.

4. Enter the current month by pressing the «**MEMORY**» button. Example: 7 x press.

5. Please follow the instructions above to set the day, hour and minutes.

6. After all settings have been made, click the «**TIME**» button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.


Note:

- If you want to change date and time, press and hold the time button down for approx. 3 seconds till Year flashes, then you can enter the new values as described above.
- 1 x press on «**MEMORY**» button or «**TIME**» button advances one operation.
- Holding the button down speeds up the procedure.

4.5. Select the User

This blood pressure monitor is designed to store 30 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.


- a) With the unit off, press and hold the «**TIME**» button for 3 seconds until the user icon in the lower left corner of the LCD screen flashes.
- b) Press the «**MEMORY**» button to toggle between users.
- c) Press the  button to make your selection.

4.6. Select Measuring Mode: Standard or Averaging Mode

This instrument enables you to select either standard (standard single measurement) or Averaging Mode (automatic triple measurement).

To select standard mode, push the switch on the side of the instrument towards «**1**», then only one measurement will be taken. If you select «**3**», the unit will inflate and deflate three times resulting in one final average measurement.

Averaging Mode Technology

- In Averaging Mode, three measurements are automatically taken in succession and the result then automatically analyzed and displayed. Because your blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.
- By pushing mode switch to 3, the selected Averaging Mode appears in the display as the  symbol.
- The bottom, right hand section of the display shows a 1, 2 or 3 to indicate which of the 3 measurements is currently being taken.
- There is a break of 15 seconds between the measurements (15 seconds are adequate according to «Blood Pressure Monitoring, 2001, 6:145-147» for oscillometric instruments). A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.
- Do not remove the cuff between measurements.

5. Carrying out a Measurement

5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try to find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Take measurements regularly at the same time of day, since blood pressure changes during the course of the day.

5.2. Common sources of error

Note: Comparable blood pressure measurements always require the same conditions!

These are normally always quiet conditions.




- All efforts by the patient to support their arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (higher) than the heart, an false reading will be obtained. Each 15 cm difference in height results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer.

Note: Only use clinically approved Cuffs!

- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.


5.3. Fitting the cuff

- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm.
- Tighten the free end of the cuff and close the cuff with the closer.
- The cuff should be snug on your upper arm, but not too tight. Any clothing which restricts the arm (e.g. pullover) must be taken off.
- Lay your arm on a table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.

Note:

If it is not possible to fit the cuff to the left arm, it can also be placed on the right. However all measurements should be made using the same arm.

**5.4. Measuring procedure**

After cuff has been appropriately positioned, the measurement can begin:

- a) Press  button. The pump begins to inflate the cuff.
Rising pressure in the cuff is shown in the display.


- b) After suitable inflation pressure has been reached, pump stops and pressure gradually falls, with cuff pressures being displayed.

In case that inflation pressure is not sufficient, the monitor automatically re-inflates to a higher level.

- c) When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.

- d) A longer beep is sounded when the measurement has been completed.
The systolic, and diastolic blood pressures and pulse rate now appear in the display.

- e) The appearance of this symbol signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements.

NOTE: We recommend contacting your physician if you see this indicator frequently.

- f) The measurement readings remain on the display until you switch off the instrument.
If no button is pressed for a period of 1 minute the instrument switches itself off in order to preserve the batteries.

- g) When the unit is set to Average Mode setting, generally 3 separate measurements will take place in succession and calculates your detected blood pressure value. There will be 15 seconds resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.


In case that the single data of each cycle differ too much from each other, a fourth measurement is performed before the result will be displayed. In rare cases the blood pressure is such unstable that even after four measurements the data vary too much. In this case «**ERR 6**» is shown and no result can be given. If one of measurement causes an error message it is repeated.

5.5. Memory-storage and recall of measurements

The blood-pressure monitor automatically stores the last 30 measurement values for each of 2 users. By pressing «**MEMORY**» button, the screen will show the most recent readings along with time and date.

Each time you press the memory button an earlier measurement will be displayed.

If an irregular heartbeat interval was detected in any measurements stored in memory, the  icon will be displayed.  indicates measurement was taken under average mode.

Memory reading


Further information

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes in a relaxed position, sitting or lying, before you repeat a measurement.

5.6. Discontinuing a measurement



If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the  button can be pressed at any time.

The device immediately lowers the cuff pressure automatically.

5.7. Delete memories

To delete stored memory readings, hold on memory button down till LCD shows a symbol of «**CL**» flashing. Then press «**Memory**» button again, you will hear 3 short beeps indicating memories have been deleted.

Before you delete all memory readings stored in, make sure you won't need refer to them at a later date. Keeping a written record is prudent and may provide additional information for your doctor's visit.

6. Error Messages/Malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed. (Example: Error no. 1)

Err
1

Error No.	Possible cause(s)
ERR 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged after the systolic blood pressure was measured. Further possible cause: No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
ERR 6	Single data differ too much during Average Mode even after 4 cycles. No average result can be displayed.
HI	The pressure in the cuff is too high (over 300mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* Please consult your doctor, if this or any other problem occurs repeatedly.

Other possible malfunctions and possible solutions – If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Disfunctions	Solution
The display remains empty when the instrument is switched on although batteries have been installed.	<ol style="list-style-type: none"> 1. Check the polarity of the batteries (+/-) 2. If the display is erratic or unusual, remove the batteries and re-install new batteries.
The pressure does not rise although the pump is running.	<ul style="list-style-type: none"> • Check the connection of the cuff tube and connect properly if necessary.
The instrument frequently fails to measure blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"> 1. Ensure the cuff fits correctly. 2. Make sure that the cuff is not fitted too tightly. Ensure tight clothing such as a rolled-up sleeve isn't exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary. 3. Measure blood pressure again in complete peace and quiet.
Every measurement produces a different value although the instrument functions normally and the values displayed are normal	<ul style="list-style-type: none"> • Please read the following information and the points listed under (Common sources of error.) Repeat the measurement.
The Blood pressure values measured differ from those measured by the doctor.	<ul style="list-style-type: none"> • Record the daily development of the values and consult your doctor about them. • Please note: Blood pressure fluctuates continually so successive measurements will show some variability.

Further information


The level of blood pressure is subject to fluctuations even in healthy people. It is important to compare measurements taken under the same conditions and at the same time of day. (Quiet conditions)!

If you have any questions regarding the use of this blood pressure monitor, please ask your dealer or pharmacist for the Microlife Service representative in your country. The Service Team will be happy to help you. **Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all guarantee claims!

7. Care and Maintenance, Recalibration

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.


- b) The cuff contains a sensitive air-tight bubble. Handle the cuff carefully and avoid all types of straining through twisting or buckling.

- c) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Carefully remove spots on the cuff with a damp cloth and soapsuds. **Do not wash the cuff in a clothes or dishwasher!**


- d) Handle the tube carefully. Avoid stress from pulling, pinching and sharp edges.

- e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

- f) **Never open the device!** This can negatively impact the calibration! (accuracy)

Periodical re-calibration

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your device by an authorized dealer every **2 years**.

Your factory authorized dealer would be pleased to provide more information regarding calibration.

8. Guarantee

This blood pressure monitor is **guaranteed for 5 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, damage from leaking batteries, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

9. Technical Specifications and reference to standards

Weights:	480 g (with batteries & cuff)
Size:	148 (I) x 126 (L) x 60 (H) mm
Storage conditions:	-20 à +55 °C; 15 à 95 % relative humidity maximum
Operating conditions:	10 à 40 °C; 15 à 95 % relative humidity maximum
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	20 à 280 mmHg
Pulse:	40 à 200 battements par minute
Cuff pressure display range:	0-299 mmHg
Memory:	2 users, each 30 memories
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse ± 5 % of the reading
Power source:	DC 6V / 600mA a) 4 dry alkaline cells (batteries), 1.5V size AA b) AC adapter 6 V DC 600 mA (optional)
Battery lifetime :	approx. 920 measurements (using new batteries)
Accessories:	M-L cuff for arm circumference 22-42 cm
Device standard:	EN 1060-1/-3/-4; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11 This device complies with the requirements of the Medical Device Directive 93/42/EEC.
IP Class:	IP20
Expected service life:	Device: 5 years or 10'000 measurements Accessories: 2 years

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S.A.V. LABORATOIRE MARQUE VERTE, A L'ADRESSE INDIQUEE CI-DESSOUS
To return with the blood pressure monitor to S.A.V. LABORATOIRE MARQUE VERTE
Mit dem Blutdruckmessgerät zurückzuschicken an S.A.V. LABORATOIRE MARQUE VERTE
Om met de bloeddrukmonitor aan S.A.V. LABORATOIRE MARQUE VERTE terug te keren